

## Membranes for Transfer and Immobilization

BioTrace™, Biodyne®, FluoroTrans® and UltraBind™ Membranes

**A variety of membrane chemistries for sensitive detection and consistent results in all applications and detection systems**

- **BioTrace NT membrane**—pure unsupported nitrocellulose is used for colony and plaque lifts and has low burn-through in protein transfer applications.
- **BioTrace PVDF membrane**—versatile membrane with broad chemical resistance, ideally suited for protein transfers.
- **Biodyne membranes**—intrinsically hydrophilic Nylon 6,6 membranes provide high sensitivity and low background for enhanced resolution. Ideal for nucleic acid blots and protein ELISA tests.
- **FluoroTrans PVDF membranes**—increased sensitivity, ideal for a wide variety of protein analysis applications including sequencing and western transfers.



- **UltraBind membrane**—affinity membrane is recommended for covalent protein binding.

### Applications

- Nucleic acid and protein transfer and detection:
  - Northern, Southern, and Western transfers
  - Colony and plaque lifts
  - Replica plating
  - Dot/slot blots
  - DNA fingerprinting
- Protein sequencing
- Solid phase ELISAs
- Affinity separations
- Macroarrays
- Microarrays

### Ordering Information

Description	Packaging	Product Numbers										
		BioTrace NT Membrane	BioTrace PVDF Membrane	Biodyne A Membrane, 0.2 µm	Biodyne A Membrane, 0.45 µm	Biodyne A Membrane, 1.2 µm	Biodyne B Membrane, 0.45 µm	Biodyne C Membrane, 0.45 µm	Biodyne Plus Membrane, 0.45 µm	FluoroTrans W Membrane, 0.2 µm	FluoroTrans Membrane, 0.2 µm	UltraBind Membrane, 0.45 µm
82 mm discs	50/pkg	66487			60102		60202	60316	60402			
85 mm discs	50/pkg	66595			60103		60203	60317	60403			
132 mm discs	50/pkg	66518			60104		60204	60318	60404			
137 mm discs	50/pkg	66488			60105		60205	60319	60405			
7 x 8.5 cm sheets	10/pkg	66593	66594		60101		60201	60315	60401		PVM020C-160	
7 x 9 cm sheets	10/pkg									BSP0158		
8.5 x 9 cm sheets	20/pkg										PVM020C-195	
10 x 15 cm sheets	10/pkg									BSP0157	PVM020C-1015	
13 x 14 cm sheets	10/pkg										PVM020C-196	
20 x 20 cm sheets	10/pkg	66489	66542		60100		60200	60314	60400	BSP0159	PVM020C-2020	66544
30 cm x 3 m roll	1/pkg	66485	66543	60113	60106	60108	60207	60320	60406			66545
20 cm x 1 m roll	1/pkg		66547				60209					
20 cm x 3 m roll	1/pkg				60120		60208					
3.3 m roll	1/pkg									BSP0161	PVM020C-099	

In addition to standard sizes, these membranes can be cut to size to suit your specifications. For information on special-sized cuts, call your local Pall Life Sciences office.

## Transfer and Affinity Membrane Selection Guide

Pall Life Sciences offers membranes for use in transfer and immobilization procedures. These membranes can be used for nucleic acid and protein applications and are compatible with radioactive, as well as nonradioactive detection systems.

Product	Biodyne® A Membrane	Biodyne B/Plus Membrane	Biodyne C Membrane
<b>Description</b>	Amphoteric Nylon 6,6	Positively-charged Nylon 6,6	Negatively-charged Nylon 6,6
<b>Works best for:</b> <b>Also suited for:</b>	Colony/Plaque Lifts, DNA and RNA Transfers  Gene Probe Assays, DNA Fingerprinting, Nucleic Acid Dot/Slot Blots, Replica Plating, ELISAs	DNA and RNA Transfers, Multiple Rebindings  DNA Fingerprinting, Nucleic Acid Dot/Slot Blots, Colony/Plaque Lifts (Biodyne B membrane), Replica Plating (Biodyne B membrane)	Reverse Dot Blots  Protein Immobilization, Affinity Purification, ELISAs
<b>Advantages</b>	<ul style="list-style-type: none"> <li>– High sensitivity</li> <li>– Low background</li> <li>– Net charge can be controlled by changing pH</li> <li>– Ability to strip and reprobe</li> </ul>	<ul style="list-style-type: none"> <li>– Positive charge over broad pH range</li> <li>– Highest sensitivity for nucleic acid applications (Biodyne B membrane)</li> <li>– Ability to strip and reprobe</li> </ul>	<ul style="list-style-type: none"> <li>– Negative charge over broad pH range</li> <li>– Surface carboxyl groups can be derivatized</li> <li>– Ability to strip and reprobe</li> </ul>
<b>Binding Interaction</b>	Hydrophobic & Electrostatic	Hydrophobic & Electrostatic	Hydrophobic & Electrostatic
<b>Method of Immobilization</b>	UV Crosslink Baking	Can be baked or UV crosslinked, although not required	Derivatization
<b>Detection Methods</b>	Radiolabeled Probes, Enzyme-antibody Conjugates <ul style="list-style-type: none"> <li>– Chemiluminescent</li> <li>– Chromogenic</li> </ul>	Radiolabeled Probes, Enzyme-antibody Conjugates <ul style="list-style-type: none"> <li>– Chemiluminescent</li> <li>– Chromogenic</li> <li>– Chemifluorescent (Biodyne Plus Membrane)</li> </ul>	Radiolabeled Probes, Enzyme-antibody Conjugates <ul style="list-style-type: none"> <li>– Chromogenic</li> </ul>

Product	BioTrace™ NT Membrane	BioTrace PVDF Membrane	FuoroTrans® Membrane	UltraBind™ Membrane
<b>Description</b>	100% Pure Nitrocellulose	Polyvinylidene Fluoride	Polyvinylidene Fluoride	Modified Polyethersulfone
<b>Works best for:</b> <b>Also suited for:</b>	Colony/Plaque Lifts  Nucleic Acid and Protein Transfers, Protein Dot/Slot Blots	Protein Transfers  Protein Dot/Slot Blots	Western Transfers (FuoroTrans W) N-terminal Protein Sequencing (FuoroTrans PVDF)	Solid-phase ELISAs  Affinity Chromatography, Hybridoma Screening
<b>Advantages</b>	<ul style="list-style-type: none"> <li>– Excellent strength</li> <li>– No support fabric</li> <li>– No detergents added</li> <li>– 100% pure nitrocellulose</li> </ul>	<ul style="list-style-type: none"> <li>– Chemical resistance</li> <li>– No discoloration</li> <li>– Nonflammable</li> <li>– High strength</li> </ul>	<ul style="list-style-type: none"> <li>– Strong protein binding</li> <li>– Sensitive detection</li> <li>– Very low burn-through</li> <li>– Good chemical compatibility</li> </ul>	<ul style="list-style-type: none"> <li>– Covalent binding</li> <li>– No preactivation required</li> <li>– High protein-binding capacity</li> </ul>
<b>Binding Interaction</b>	Hydrophobic & Electrostatic	Hydrophobic	Hydrophobic	Covalent
<b>Method of Immobilization</b>	UV Crosslink Baking (Vacuum Oven)			Direct Spotting Perfusion
<b>Detection Methods</b>	Radiolabeled Probes, Direct Stain, Fluorescence, Enzyme-antibody Conjugates <ul style="list-style-type: none"> <li>– Chemiluminescent</li> <li>– Chromogenic</li> </ul>	Direct Stain, Enzyme-antibody Conjugates <ul style="list-style-type: none"> <li>– Chemiluminescent</li> <li>– Chromogenic</li> </ul>	Direct Stain with Coomassie blue, Amido black, Ponceau S, and colloidal gold (FuoroTrans W membrane). Enzyme-antibody Conjugates <ul style="list-style-type: none"> <li>– Chemiluminescent</li> <li>– Chromogenic</li> </ul>	Radiolabeled Probes, Enzyme-antibody Conjugates <ul style="list-style-type: none"> <li>– Chromogenic</li> </ul>

## Biodyne® Transfer Membranes

- High sensitivity and low background for enhanced detection and resolution.
- Do not crack, shrink, or tear when subjected to multiple cycles of hybridization, stripping, and reprobing.
- Intrinsically hydrophilic for easy wetting.
- Superior performance with radioactive (Biodyne B membrane) and nonradioactive (Biodyne A membrane) detection systems.
- Ideal for nucleic acid detection.

### Applications

Four chemistries provide versatile adsorption properties:

1. **Biodyne A Membrane: Amphoteric Nylon 6,6.**  
Membrane zeta potential can be modulated by changes in pH. Ideal for single probe or multiple rehybridizations, and applications where background is troublesome.
2. **Biodyne B Membrane: Positively-charged Nylon 6,6.**  
Pore surfaces are populated by a high density of quaternary ammonium groups. Our highest sensitivity nylon membrane for nucleic acid applications.
3. **Biodyne C Membrane: Negatively-charged Nylon 6,6.**  
Can be derivatized by coupling reactions through the carboxyl groups on the pore surfaces.
4. **Biodyne Plus Membrane: Positively-charged Nylon 6,6 with an extremely high isoelectric point.**  
With certain nonradioactive detection systems, it is more sensitive than Biodyne A membrane while exhibiting lower background than Biodyne B membrane.

### Specifications

#### Media

Nylon 6,6

#### Typical Thickness

6.0 mils ± 0.5 mils

#### Pore Sizes


0.2, 0.45, or 1.2 µm

#### Solvent Compatibility


Resistant to common solvents such as acetone, alcohol, chlorinated aliphatic hydrocarbons, formamide, 2 M NaOH, DMSO, and dimethylformamide. Not compatible with concentrated formic acid (> 50%), HCl (> 4 M), oxidizing agents, and long exposures (days to weeks) at pH < 2.

### Performance

Biodyne B Membrane Withstands Multiple Cycles of Stripping and Reprobing



Superior Fluorescent Detection of DNA Using Biodyne Plus Membrane



## BioTrace™ PVDF Transfer Membrane

- Versatile membrane for nucleic acid and protein transfers.
- Broad compatibility with commonly-used solvents.
- Low background with chemiluminescent detection systems.

---

### Applications

- Western transfers
- Southern transfers

---

### Specifications

#### Media

Polyvinylidene fluoride

#### Typical Thickness

147  $\mu\text{m}$  (5.8 mils)

#### Pore Size

0.45  $\mu\text{m}$

#### Tensile Strength

28 bar (410 psi)


#### Solvent Compatibility

Resistant to methanol, phenol, and chloroform. Also resistant to 10% dimethyl sulfoxide, 15% acetic acid, 70% formic acid, 25% triethylamine, 1 N NaOH, and 1 N KOH.

---

### Performance

Western Transfer to BioTrace PVDF Membrane



## BioTrace™ NT Transfer Membrane

- Pure unsupported nitrocellulose membrane is ideal for colony/plaque lifts and protein transfers.
- Strong and durable, less likely to tear or crack than competitor nitrocellulose.
- High binding capacity for proteins and nucleic acids.
- Lower protein burn-through than competitors in electrophoretic transfers.

---

### Applications

- Colony/plaque lifts
- Protein transfers

---

### Specifications

#### Media

Nitrocellulose

#### Typical Thickness

145  $\mu\text{m}$  (5.7 mils)

#### Pore Size

0.2  $\mu\text{m}$


#### Protein Binding Capacity

209  $\mu\text{g}/\text{cm}^2$

---

### Performance

BioTrace NT Membranes Exhibit Low Protein Burn-through



## FluoroTrans® PVDF Membrane

- Sensitive protein detection with low background and very low burn-through.
- Membranes provide high surface area for strong hydrophobic interactions and typically adsorb 50% more protein than nylon or nitrocellulose.
- FluoroTrans W membrane is optimized for Western transfer applications.
- FluoroTrans PVDF membrane is optimized for N-terminal protein sequencing.

---

### Applications

FluoroTrans W Membrane:

- Western transfers
- Southern transfers

FluoroTrans Membrane:

- N-terminal protein sequencing

- FluoroTrans media have high tensile strength and will not tear, crack, or curl during handling. This allows for easy removal of target bands for protein sequencing applications.

---

### Specifications

#### Media

Hydrophobic polyvinylidene fluoride

#### Pore Size

0.2 µm

#### Chemical Compatibility

Resistant to acetone, DMSO, dimethyl formamide, methanol, trifluoroacetic acid, and triethylamine.

---

### Performance

FluoroTrans Membrane has Excellent Sensitivity, Signal, and Background in Western Transfers



# UltraBind™ Affinity Membrane

- Modified polyethersulfone (PES) membrane for covalent protein binding.
- Proteins can be efficiently attached without prior membrane derivitization.

## Applications

- ELISA
- Affinity separation

## Specifications

### Media

Modified polyethersulfone with aldehyde surface chemistry

### Pore Size

0.45  $\mu\text{m}$

### Typical Thickness


152  $\mu\text{m}$  (6 mils)

### Typical IgG Binding Capacity


135  $\mu\text{g}/\text{cm}^2$

## Performance

UltraBind Membrane Binds and Retains Proteins



Antigen Detection (dot blot ELISA) with UltraBind Membrane



## Complementary Products

- **Centrifugal Devices** provide precise, rapid processing of the following sample volumes:

Device	Sample Volumes
Nanosep <sup>®</sup> Device	50 to 500 µL
Microsep <sup>™</sup> Device	500 µL to 3.5 mL
Macrosep <sup>®</sup> Device	1 mL to 15 mL
Jumbosep <sup>™</sup> Device	15 mL to 60 mL

- **AcroWell<sup>™</sup> 96- and 384-well Filter Plates** with BioTrace Membranes exhibit high binding capacities for proteins and nucleic acids.
- **AcroPrep<sup>™</sup> 96- and 384-well Filter Plates** can be used for a variety of molecular biology, combinatorial chemistry, and screening applications.
- **Vivid<sup>™</sup> Gene Array Slides** feature a unique membrane construction that allows high signal-to-noise ratios, requires less template, and provides consistent results. Protocols are easy to follow with simple immobilization steps.

## Technical Literature

- Discover Endless Potential: Products for Genomics, Proteomics, and Drug Discovery Brochure, Pall Life Sciences, PN33286
- Transfer and Detection Procedures for Pall Life Sciences Membranes and Kits, Pall Life Sciences, PN33167
- Explore the Possibilities: High Throughput Separation, Purification, and Detection Technologies Brochure, Pall Life Sciences, PN33252



Life Sciences

**Pall Life Sciences**  
600 South Wagner Road  
Ann Arbor, MI 48103-9019 USA

800.521.1520 toll free in USA  
734.665.0651 phone  
734.913.6114 fax

### Australia – Lane Cove, NSW

Tel: 02 9428-2333  
1800 635-082 (in Australia)  
Fax: 02 9428-5610

### Austria – Wien

Tel: 043-1-49 192-0  
Fax: 0043-1-49 192-400

### Canada – Ontario

Tel: 905-542-0330  
800-263-5910 (in Canada)  
Fax: 905-542-0331

### Canada – Québec

Tel: 514-332-7255  
800-435-6268 (in Canada)  
Fax: 514-332-0996

### 800-808-6268 (in Canada)

**China – P. R., Beijing**  
Tel: 86-10-8458 4010  
Fax: 86-10-8458 4001

### France – St. Germain-en-Laye

Tel: 01 30 61 39 92  
Fax: 01 30 61 58 01  
Lab-FR@pall.com

### Germany – Dreieich

Tel: 06103-307 333  
Fax: 06103-307 399  
Lab-DE@pall.com

### India – Mumbai

Tel: 91-22-5956050  
Fax: 91-22-5956051

### Italy – Milano

Tel: 02-47796-1  
Fax: 02-47796-394  
or 02-41-22-985

### Japan – Tokyo

Tel: 3-3495-8319  
Fax: 3-3495-5397

### Korea – Seoul

Tel: 2-569-9161  
Fax: 2-569-9092

### Poland – Warszawa

Tel/Fax: 22-835 83 83

### Russia – Moscow

Tel: 095 787-76-14  
Fax: 095 787-76-15

### Singapore

Tel: (65) 389-6500  
Fax: (65) 389-6501

### Spain – Madrid

Tel: 91-657-9876  
Fax: 91-657-9836

### Sweden – Lund

Tel: +46 (0)46 158400  
Fax: +46 (0)46 320781

### Switzerland – Basel

Tel: 061-638 39 00  
Fax: 061-638 39 40

### Taiwan – Taipei

Tel: 2-2545-5991  
Fax: 2-2545-5990

### United Kingdom – Portsmouth

Tel: 023 92 302600  
Fax: 023 92 302601

Lab-UK@pall.com

Visit us on the Web at [www.pall.com/lab](http://www.pall.com/lab)

E-mail us at [Lab@pall.com](mailto:Lab@pall.com)

© 2003, Pall Corporation. Pall, , AcroPrep, AcroWell, Biodyne, BioTrace, FluoroTrans, Jumbosep, Macrosep, Microsep, Nanosep, UltraBind and Vivid are trademarks of Pall Corporation. ® indicates a trademark registered in the USA. **Filtration. Separation. Solution.<sup>sm</sup>** is a service mark of Pall Corporation. FluorImager is a registered trademark of Molecular Dynamics. PlateMate is a trademark of Matrix Corporation. LumiGLO is a registered trademark of Kirkegaard & Perry Laboratories, Inc.

**Filtration. Separation. Solution.<sup>sm</sup>**